



**HERVIK**<sup>®</sup>  
Healthcare Pvt Ltd



**Broucher:**  
**M/s. Hervik Healthcare Private Limited.**

e-mail: [hervikhealthcare@gmail.com](mailto:hervikhealthcare@gmail.com)  
website: [www.hervikhealthcare.com](http://www.hervikhealthcare.com)  
Customer care: +91 9985 979 770.

# Disposable Blood Transfusion / administration Sets

(Ribbon, Blister & Pouch Packed “Higgs & Hervik” brands Disposable Blood Transfusion / administration Sets).

Disposable Blood Transfusion set it is used to administer blood from a container to a patient's vascular system through a needle or catheter inserted into a vein.



#### PRODUCT SPECIFICATIONS:

- Specially designed for Intravenous Infusion of Blood and Blood derivatives.
- Drip formation part designed as per IS Standards for the distilled water flow rate 20 drops = {1.0 + 1 ml (1.0 + 0.1 g).
- Transparent Drip Chamber easy priming and easy monitoring of flow rate.
- 10cm<sup>2</sup> area Blood filter of 200 micron.
- 304 Imported cannula tube for needle & it's siliconized with imported silicon.
- Medical grade, Nontoxic PVC Tubing for kink free effect.
- Sterile, disposable, non-pyrogenic, individually packed.
- BLOOD TRANSFUSION SETS CONFORMS AS PER IS 9824 (PART - 31996).
- Individually pack in each unit. In both (poly pack & Ribbon pack).

#### PACKING OPTIONS:

- Individually packing: Medical Grade paper pouch with standard plastic Film/Ribbon Pack
- Inner packing: Corrugated carton box
- Outer packing: Corrugated carton box

#### PRECAUTION:

- Before use, check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.
- Infusion Set is designed for single use only
- Measured Set is designed for single use only
- It is recommended that set is changed every 24 hours

#### WARNING:

- Do not attempt to re-sterilize and reuse.
- Not to be used for "Blood and Blood derivatives"

#### PRODUCT RANGE:

- With Air vent and without air vent and other variants are available.



**Our Prestigious Brand:**



**M/s. Hervik Healthcare Private Limited**

Manufacturers of Disposable Medical Devices.

Reg. Office: Door Number. 8-202, Jamandla Street,  
Chandragiri, Chittoor (Dist),  
Andhra Pradesh - 517101, India.

Factory Address: Survey No: 600/B,  
NH-140, Pichinaidupalli Road,  
Chandragiri, Chittoor (Dist),  
Andhra Pradesh - 517101, India.

Contact Mobile Number: +91 9985 979 770  
Email: [hervikhealthcare@gmail.com](mailto:hervikhealthcare@gmail.com)  
Website: [www.hervikhealthcare.com](http://www.hervikhealthcare.com)