



**HERVIK**<sup>®</sup>  
Healthcare Pvt Ltd



**Broucher:**  
**M/s. Hervik Healthcare Private Limited.**


e-mail: [hervikhealthcare@gmail.com](mailto:hervikhealthcare@gmail.com)  
website: [www.hervikhealthcare.com](http://www.hervikhealthcare.com)  
Customer care: +91 9985 979 770.


# Hypodermic Needles - [DISPOSABLE NEEDLES]



(Blister Packed "Higgs & Hervik" brands Hypodermic Disposable needles).

Hypodermic Disposable Needles are manufactured using premium quality, stainless steel with the latest Japanese Technology, "Higgs" needles are ground on high-precision machines and assembled and packed on totally automated machines with vision camera quality controls for enabling superior sharpness to ensure (smooth penetration) injections with minimal pain.





**PRODUCT SPECIFICATIONS:**

- Disposable Needles are made from imported 304 stainless steel Cannula
- Available as three-face-bevel & Ultra-sharp lancet point
- Siliconized, minimizing drag during penetration and withdrawal.
- Needle point is checked by vision camera system.
- Non-toxic and Non-pyrogenic
- Disposable needles are sterilized with Ethylene Oxide (EO) gas.
- Manufactured as per IS : 10654 /ISO:7864 / ISO:6009

**PRODUCT PACKING:**

- Individually packing: The disposable needles are packed in blister packs having easily peelable, fibre-free, imported, medical-grade paper to ensure reliability of sealing and product sterility.
- Inner packing: Corrugated carton box
- Outer packing: Corrugated carton box

**PRECAUTION:**

- Before use, check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.

**WARNING:**

- Do not attempt to re-sterilize and reuse.

**PRODUCT RANGE:**

S. No	Needle Gauge	Needle Hub Color	Needle Length
1	16G	White	As per customer requirement, we will provide any gauge and any length (0.5", 1", 1.25" & 1.5") of needles.
2	18G	Pink	
3	21G	Green	
4	22G	Black	
5	23G	Blue	
6	24G	Medium Purple	
7	26G	Brown	



**Our Prestigious Brand:**



**M/s. Hervik Healthcare Private Limited**

Manufacturers of Disposable Medical Devices.

Reg. Office: Door Number. 8-202, Jamandla Street,  
Chandragiri, Chittoor (Dist),  
Andhra Pradesh - 517101, India.

Factory Address: Survey No: 600/B,  
NH-140, Pichinaidupalli Road,  
Chandragiri, Chittoor (Dist),  
Andhra Pradesh - 517101, India.

Contact Mobile Number: +91 9985 979 770  
Email: [hervikhealthcare@gmail.com](mailto:hervikhealthcare@gmail.com)  
Website: [www.hervikhealthcare.com](http://www.hervikhealthcare.com)