



HERVIK[®]
Healthcare Pvt Ltd



Broucher:
M/s. Hervik Healthcare Private Limited.

e-mail: hervikhealthcare@gmail.com
website: www.hervikhealthcare.com
Customer care: +91 9985 979 770.

"Why HERVIK...!"

We Strive to be the Most Trusted medical devices Manufacturing Solutions Provider.



CREATIVITY



- **SOLUTIONS:** Our persistent drive to tailor SOLUTIONS for each customer success. That we strongly believe builds lasting partnerships with HERVIK.
- **INGENUITY:** We are innovators at heart; and we think BIG.
- **INTELLIGENCE:** This applied data INTELLIGENCE is also what drives our ability to help customers win in an increasingly medical devices and healthcare sectors.
- **RELIABILITY:** "Do the Right Thing" has been our mantra from the beginning. Our partnerships, experience, business model and RELIABILITY make us the trusted pair.
- **IMPACT:** We are a trusted brand that exists to serve our customers. But, as a responsible global citizen, we are proud to champion a common purpose to make a positive IMPACT for our healthcare sector.
- **OPPORTUNITY:** Everyone's voice is heard, and we work together to bring OPPORTUNITY to each other, our customers and our shareholders.

About Hervik Healthcare

What makes us and Who we are today...



Hervik Healthcare Private Limited (HHPL) was founded in 2020 and created into manufacturing and marketing (globally) of medical devices and disposables. HHPL markets its products under “HERVIK & HIGGS” brands and we are one of the finest producers of syringes in India with a great installed capacity that aims to maintain International Standards at all costs; thereby offering cutting-edge medical consumables to a highly demanding global market. Driven by a dedicated and meticulous production team.

Hervik Healthcare Private Limited (HHPL) is having fully equipped manufacturing plant in Chandragiri, Chittoor (Dist.) of Andhra Pradesh, India.

The following is the current list of HHPL products are:

- Disposable Hypodermic syringes (Single use) (with all variants).
- Disposable Hypodermic syringes (Auto-disable / Safety) (Single use).
- Disposable Hypodermic Needles.
- I.V. Cannulas.
- Infusion / Transfusion sets.
- And Other medical / disposable items.

Above all products are manufacturing and marketed through an established national distribution network, to ensure ready availability even in remote areas, and effective customer service support.

Hervik Healthcare Private Limited (HHPL) marketing philosophy is to manufacture products of international standards and quality with affordable prices. Hervik Healthcare Private Limited (HHPL) manufacturing philosophy is to achieve the best quality, through a dedicated production team, in a harmonious industrial atmosphere as per requirement of ISO 13485.

Hervik Healthcare Private Limited (HHPL) it's vision and long-term goal of becoming 'Top 3' in each of its business segments & building a unique competitive advantage in India.

Manufacturing Facility

What makes us and Who we are today...





Hervik - Field of Manufacturing



Hervik Healthcare Private Limited (HHPL) is the most promising manufacturing organization that serves the healthcare industry with comprehensive services from medical devices development through commercial manufacturing. Hervik Healthcare Private Limited (HHPL) has fully-integrated plant in the town of Chandragiri, Tirupati, Andhra Pradesh, India. This plant is designed to manufacture the entire range of Hervik Healthcare Private Limited (HHPL) products. This plant is also equipped with state-of-art design and engineering facilities for developing new products and innovative variants of existing products, such as Disposable Syringes, Needles, (AD) Syringes, I.V. Cannulas and other infusion sets. Our Hervik Healthcare Plant run 24x7.

Hervik Healthcare Private Limited (HHPL) is Operated by professionals, HH is a known for being an employee focused Company, Hervik Healthcare (HH) places utmost focus on personal development of human skills, marketing techniques, and safety standards for its work force. The organization is dedicated to its moral obligations towards the environment and community as a whole.





Hypodermic Disposable syringes

(Ribbon and Blister Packed “Higgs & Hervik” brands Hypodermic Disposable Syringes with needle and without needle).



3part Disposable Syringes (Luer Mount/slip)



2 Part Disposable Syringes (Luer Mount/slip)



3 Part Disposable Syringes (Luer Lock)



PRODUCT SPECIFICATIONS:

- Luer Slip or mount / luer lock Single use syringes with needle & without needle.
- Manufactured from medical grade virgin polypropylene.
- Prominent graduation ensures dosage accuracy and avoids wastage of valuable medicine.
- Latex free & thermoplastic Elastomer gasket is inert that provides minimum friction during movement and prevents leakage and back flow.
- Sharp triple faceted bevelled tip needle for minimum pain and extremely smooth penetration.
- Sterile, Disposable (Single use only), Non-Pyrogenic, Latex Free & Individually packed.
- Manufactured as per IS: 10258 / ISO: 7886-1.

PRODUCT PACKING:

- Individually packing: Ribbon Pack in tough laminated film, highly resistant to puncture and tearing.
- Individually packing: Blister Pack in medical grade paper and film making it safe and suitable for long-term storage.
- Inner packing: Corrugated carton box
- Outer packing: Corrugated carton box
- For more packaging details please visit www.hervikhealthcare.com.

PRECAUTION:

- Before use, check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.

WARNING:

- Do not attempt to re-sterilize and reuse.

PRODUCT RANGE:

Capacity	Types of syringes	Model of syringe	Standard Needle Dimensions
0.5 ml	3 Part Syringes	luer lock / luer slip	Needle: 24G, 26G, 28G, 30G, 31G
1 ml	3 Part Syringes	luer lock / luer slip	Needle: 24G, 26G, 28G, 30G, 31G
2 ml	Both 2 part & 3 Part Syringes	luer lock / luer slip	Needle: 18G, 20G, 22G, 23G, 24G
3 ml	3 Part Syringes	luer lock / luer slip	Needle: 22G, 23G, 24G
5 ml	Both 2 part & 3 Part Syringes	luer lock / luer slip	Needle: 18G, 20G, 21G, 22G, 23G, 24G
10 ml	Both 2 part & 3 Part Syringes	luer slip	Needle: 16G, 18G, 20G, 21G, 22G, 23G
20 ml	3 Part Syringes	luer slip	Needle: 18G, 20G, 21G

Advantages:

- Note: As per customer requirement, we will provide any gauge and any length (0.5", 1", 1.25" & 1.5") of needles.


Hypodermic Auto-disable (AD)/ Safety Disposable syringes



(Ribbon and Blister Packed "Higgs & Hervik" brands Hypodermic Auto-disable / safety Disposable Syringes with needle and without needle).



Auto-disable (AD) syringes automatically become disabled after one use. The "Higgs" syringe is a single-use, disposable, auto-disable syringe with a safety plunger that breaks off after a single use. The "Higgs" syringes can be purchased with a detachable needle. Auto disable syringes are meets ISO 7886-3 & ISO7886-4 Standards for A.D. Syringes.



PRODUCT SPECIFICATIONS:

- Auto-disable (AD) syringes automatically become disabled after one use.
- AD Syringes are available with needle & without needle.
- Manufactured from medical grade virgin polypropylene.
- Prominent graduation ensures dosage accuracy and avoids wastage of valuable medicine.
- Latex free & thermoplastic Elastomer gasket is inert that provides minimum friction during movement and prevents leakage and back flow.
- Sharp triple faceted bevelled tip needle for minimum pain and extremely smooth penetration.
- Sterile, Disposable (Single use only), Non-Pyrogenic, Latex Free & Individually packed.
- Manufactured as per IS: 10258 / ISO: 7886-1.

PRODUCT PACKING:

- Individually packing: Ribbon Pack in tough laminated film, highly resistant to puncture and tearing.
- Individually packing: Blister Pack in medical grade paper and film making it safe and suitable for long-term storage.
- Inner packing: Corrugated carton box
- Outer packing: Corrugated carton box
- For more packaging details please download e-catalogue.

LOCK (DISABLE SYSTEM):

- Presence of loop & lock mechanism which automatically breaks the plunger, if reuse is attempted.

PRECAUTION:

- Before use, check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.
- Do not push the plunger forward before withdrawing medicament from the vial otherwise plunger will get locked.

WARNING:

- Do not attempt to re-sterilize and reuse.

PRODUCT RANGE:

Capacity	Types of syringes	Model of syringe	Standard Needle Dimensions
0.5 ml	3 Part Syringes	luer slip / luer lock	Needle: 24G, 26G, 28G, 30G, 31G
1 ml	3 Part Syringes	luer slip / luer lock	Needle: 24G, 26G, 28G, 30G, 31G
2 ml	3 Part Syringes	luer slip / luer lock	Needle: 18G, 20G, 22G, 23G, 24G
3 ml	3 Part Syringes	luer slip / luer lock	Needle: 22G, 23G, 24G
5 ml	3 Part Syringes	luer slip / luer lock	Needle: 18G, 20G, 21G, 22G, 23G, 24G
10 ml	3 Part Syringes	luer slip	Needle: 16G, 18G, 20G, 21G, 22G, 23G

Advantages:

- Note: As per customer requirement, we will provide any gauge and any length (0.5", 1", 1.25" & 1.5") of needles.

Hypodermic Disposable Insulin / Tuberculin syringes



(Ribbon and Blister Packed “Higgs & Hervik” brands Hypodermic Disposable Insulin / Tuberculin Syringes with needle and without needle).

Our Insulin Syringes are uniquely purposed syringes. The needle is short and thin and is covered with a premium silicone layer, thus allowing it to pass through the skin fairly easily.





PRODUCT SPECIFICATIONS:

- Our Insulin syringes are siliconized ultra-fine Needle, that great quality causing Minimal Pain.
- The presence of low dead spaces thus avoiding wastage of medication and ensures accurate dosage.
- Available in both U-40 & U-100 graduated scales.
- Automated Manufacturing process (Checked by Vision Inspection Cameras)
- Cannula which is using for insulin syringes those Cannula are 100% checked

PRODUCT PACKING:

- Individually packing: Ribbon Pack in tough laminated film, highly resistant to puncture and tearing.
- Individually packing: Blister Pack in medical grade paper and film making it safe and suitable for long-term storage.
- Inner packing: Corrugated carton box
- Outer packing: Corrugated carton box
- For more packaging details please visit www.hervikhealthcare.com.

PRECAUTION:

- Before use, check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.

WARNING:

- Do not attempt to re-sterilize and reuse.

PRODUCT RANGE:

Capacity	Types of syringes	Model of syringe	Standard Needle Dimensions
0.5 ml	3 Part Syringes	luer lock / luer slip	Needle: 26G, 27G, 28G, 29G, 30G, 31G
1 ml (U40, U80 & U100 / 1ml)	3 Part Syringes	luer lock / luer slip	

Advantages:


- Note: As per customer requirement, we will provide any gauge and any length (0.5", 1", 1.25" & 1.5") of needles.


Hypodermic Needles - [DISPOSABLE NEEDLES]



(Blister Packed "Higgs & Hervik" brands Hypodermic Disposable needles).

Hypodermic Disposable Needles are manufactured using premium quality, stainless steel with the latest Japanese Technology, "Higgs" needles are ground on high-precision machines and assembled and packed on totally automated machines with vision camera quality controls for enabling superior sharpness to ensure (smooth penetration) injections with minimal pain.





PRODUCT SPECIFICATIONS:

- Disposable Needles are made from imported 304 stainless steel Cannula
- Available as three-face-bevel & Ultra-sharp lancet point
- Siliconized, minimizing drag during penetration and withdrawal.
- Needle point is checked by vision camera system.
- Non-toxic and Non-pyrogenic
- Disposable needles are sterilized with Ethylene Oxide (EO) gas.
- Manufactured as per IS : 10654 / ISO:7864 / ISO:6009

PRODUCT PACKING:

- Individually packing: The disposable needles are packed in blister packs having easily peelable, fibre-free, imported, medical-grade paper to ensure reliability of sealing and product sterility.
- Inner packing: Corrugated carton box
- Outer packing: Corrugated carton box

PRECAUTION:

- Before use, check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.

WARNING:

- Do not attempt to re-sterilize and reuse.

PRODUCT RANGE:

S. No	Needle Gauge	Needle Hub Color	Needle Length
1	16G	White	As per customer requirement, we will provide any gauge and any length (0.5", 1", 1.25" & 1.5") of needles.
2	18G	Pink	
3	21G	Green	
4	22G	Black	
5	23G	Blue	
6	24G	Medium Purple	
7	26G	Brown	

Disposable Infusion sets / Perfusion Sets

(Ribbon, Blister & Pouch Packed “Higgs & Hervik” brands Disposable Infusion sets / Perfusion Sets).



Disposable Infusion sets / Perfusion Sets are intended for use in the administration of fluids from a container into the patient's vascular system through a vascular access device.



PRODUCT SPECIFICATIONS:

- Conventional Infusion Set.
- Specially designed cylindrical collapsible drip chamber to visualize the flow rate.
- Disc type fluid filter at the bottom of chamber filters any particulate matter in the I.V fluid
- Drip formation rate is 20 drops = 1.0 + 0.1 ml (1.0 + 0.1 grams)
- 150 cm long super smooth kink resistant tubing.
- Superior quality latex bulb for additional medication.
- Micro-pin provided for Micro drip (for Micro Non-Vented & Vented sets).
- Sterile, disposable, non-pyrogenic, individually packed.
- Infusion set as per ISO 8536 - 4 : 1987 / IS 12655 (PART -4) 1988

PACKING OPTIONS:

- Individually packing: Medical Grade paper pouch with standard plastic Film/Ribbon Pack
- Inner packing: Corrugated carton box
- Outer packing: Corrugated carton box

PRECAUTION:

- Before use, check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.
- Infusion Set is designed for single use only
- Measured Set is designed for single use only
- It is recommended that set is changed every 24 hours

WARNING:

- Do not attempt to re-sterilize and reuse.
- Not to be used for "Blood and Blood derivatives"

PRODUCT RANGE:

- With Air vent and without air vent and other variants are available.

Disposable Blood Transfusion / administration Sets



(Ribbon, Blister & Pouch Packed “Higgs & Hervik” brands Disposable Blood Transfusion / administration Sets).

Disposable Blood Transfusion set it is used to administer blood from a container to a patient's vascular system through a needle or catheter inserted into a vein.



PRODUCT SPECIFICATIONS:

- Specially designed for Intravenous Infusion of Blood and Blood derivatives.
- Drip formation part designed as per IS Standards for the distilled water flow rate 20 drops = {1.0 + 1 ml (1.0 + 0.1 g).
- Transparent Drip Chamber easy priming and easy monitoring of flow rate.
- 10cm² area Blood filter of 200 micron.
- 304 Imported cannula tube for needle & it's siliconized with imported silicon.
- Medical grade, Nontoxic PVC Tubing for kink free effect.
- Sterile, disposable, non-pyrogenic, individually packed.
- BLOOD TRANSFUSION SETS CONFORMS AS PER IS 9824 (PART - 31996).
- Individually pack in each unit. In both (poly pack & Ribbon pack).

PACKING OPTIONS:

- Individually packing: Medical Grade paper pouch with standard plastic Film/Ribbon Pack
- Inner packing: Corrugated carton box
- Outer packing: Corrugated carton box

PRECAUTION:

- Before use, check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.
- Infusion Set is designed for single use only
- Measured Set is designed for single use only
- It is recommended that set is changed every 24 hours

WARNING:

- Do not attempt to re-sterilize and reuse.
- Not to be used for "Blood and Blood derivatives"

PRODUCT RANGE:

- With Air vent and without air vent and other variants are available.

Disposable I.V. CANNULAS

(Blister Packed “Higgs & Hervik” brands Disposable I.V. Cannulas).



(Blister Packed “Higgs & Hervik” brands Disposable I.V. Cannulas)



PRODUCT SPECIFICATIONS:

- Premium I.V. Cannula with wings (1 way & 2 way)
- It is provided with low pressure injection valve to prevent its back flow.
- Medicine can be injected by a syringe without needle. Cost of needle is also saved.
- Ergonomically designed: The valve is protected with easy to open Injection Port Cap.
- Colored body / Injection Port Cap for size identification.
- Optimum trim distance for pain free penetration.
- Healthy, evenly oriented, sharp, triple faceted bevelled needle ensures least trauma and reliable insertion.
- Wings specially designed for easy gripping and safe clamping. It's offer easy fixation and prevents shifting & rolling of Cannula in the patient's body.
- Stainless steel (304) needle with diameter and length adhering to all ISO standards
- Single use, latex free, sterile, non-pyrogenic, nontoxic.

PACKING OPTIONS:

- Individually packing: Medical Grade paper individually packed in blister pack
- Inner packing: Corrugated carton box
- Outer packing: Corrugated carton box

PRECAUTION:

- Before use, check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.

WARNING:

- Do not attempt to re-sterilize and reuse.



Our Prestigious Brand:



M/s. Hervik Healthcare Private Limited

Manufacturers of Disposable Medical Devices.

Reg. Office: Door Number. 8-202, Jamandla Street,
Chandragiri, Chittoor (Dist),
Andhra Pradesh - 517101, India.

Factory Address: Survey No: 600/B,
NH-140, Pichinaidupalli Road,
Chandragiri, Chittoor (Dist),
Andhra Pradesh - 517101, India.

Contact Mobile Number: +91 9985 979 770
Email: hervikhealthcare@gmail.com
Website: www.hervikhealthcare.com