



**HERVIK**<sup>®</sup>  
Healthcare Pvt Ltd



**Broucher:**  
**M/s. Hervik Healthcare Private Limited.**


e-mail: [hervikhealthcare@gmail.com](mailto:hervikhealthcare@gmail.com)  
website: [www.hervikhealthcare.com](http://www.hervikhealthcare.com)  
Customer care: +91 9985 979 770.

# Disposable I.V. CANNULAS

(Blister Packed “Higgs & Hervik” brands Disposable I.V. Cannulas).



(Blister Packed “Higgs & Hervik” brands Disposable I.V. Cannulas)



## PRODUCT SPECIFICATIONS:

- Premium I.V. Cannula with wings (1 way & 2 way)
- It is provided with low pressure injection valve to prevent its back flow.
- Medicine can be injected by a syringe without needle. Cost of needle is also saved.
- Ergonomically designed: The valve is protected with easy to open Injection Port Cap.
- Colored body / Injection Port Cap for size identification.
- Optimum trim distance for pain free penetration.
- Healthy, evenly oriented, sharp, triple faceted bevelled needle ensures least trauma and reliable insertion.
- Wings specially designed for easy gripping and safe clamping. It's offer easy fixation and prevents shifting & rolling of Cannula in the patient's body.
- Stainless steel (304) needle with diameter and length adhering to all ISO standards
- Single use, latex free, sterile, non-pyrogenic, nontoxic.

## PACKING OPTIONS:

- Individually packing: Medical Grade paper individually packed in blister pack
- Inner packing: Corrugated carton box
- Outer packing: Corrugated carton box

## PRECAUTION:

- Before use, check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.

## WARNING:

- Do not attempt to re-sterilize and reuse.



**Our Prestigious Brand:**



**M/s. Hervik Healthcare Private Limited**

Manufacturers of Disposable Medical Devices.

Reg. Office: Door Number. 8-202, Jamandla Street,  
Chandragiri, Chittoor (Dist),  
Andhra Pradesh - 517101, India.

Factory Address: Survey No: 600/B,  
NH-140, Pichinaidupalli Road,  
Chandragiri, Chittoor (Dist),  
Andhra Pradesh - 517101, India.

Contact Mobile Number: +91 9985 979 770  
Email: [hervikhealthcare@gmail.com](mailto:hervikhealthcare@gmail.com)  
Website: [www.hervikhealthcare.com](http://www.hervikhealthcare.com)