



HERVIK[®]
Healthcare Pvt Ltd



Broucher:
M/s. Hervik Healthcare Private Limited.

e-mail: hervikhealthcare@gmail.com
website: www.hervikhealthcare.com
Customer care: +91 9985 979 770.

Disposable Infusion sets / Perfusion Sets

(Ribbon, Blister & Pouch Packed “Higgs & Hervik” brands Disposable Infusion sets / Perfusion Sets).



Disposable Infusion sets / Perfusion Sets are intended for use in the administration of fluids from a container into the patient's vascular system through a vascular access device.



PRODUCT SPECIFICATIONS:

- Conventional Infusion Set.
- Specially designed cylindrical collapsible drip chamber to visualize the flow rate.
- Disc type fluid filter at the bottom of chamber filters any particulate matter in the I.V fluid
- Drip formation rate is 20 drops = 1.0 + 0.1 ml (1.0 + 0.1 grams)
- 150 cm long super smooth kink resistant tubing.
- Superior quality latex bulb for additional medication.
- Micro-pin provided for Micro drip (for Micro Non-Vented & Vented sets).
- Sterile, disposable, non-pyrogenic, individually packed.
- Infusion set as per ISO 8536 - 4 : 1987 / IS 12655 (PART -4) 1988

PACKING OPTIONS:

- Individually packing: Medical Grade paper pouch with standard plastic Film/Ribbon Pack
- Inner packing: Corrugated carton box
- Outer packing: Corrugated carton box

PRECAUTION:

- Before use, check the integrity of the sealing of the pack, if the pack is previously opened or damaged then do not use.
- Infusion Set is designed for single use only
- Measured Set is designed for single use only
- It is recommended that set is changed every 24 hours

WARNING:

- Do not attempt to re-sterilize and reuse.
- Not to be used for "Blood and Blood derivatives"

PRODUCT RANGE:

- With Air vent and without air vent and other variants are available.



Our Prestigious Brand:



M/s. Hervik Healthcare Private Limited

Manufacturers of Disposable Medical Devices.

Reg. Office: Door Number. 8-202, Jamandla Street,
Chandragiri, Chittoor (Dist),
Andhra Pradesh - 517101, India.

Factory Address: Survey No: 600/B,
NH-140, Pichinaidupalli Road,
Chandragiri, Chittoor (Dist),
Andhra Pradesh - 517101, India.

Contact Mobile Number: +91 9985 979 770
Email: hervikhealthcare@gmail.com
Website: www.hervikhealthcare.com